

Product: Tissue Grossing Table with Downdraft Extraction, no sink
Cat.-No.: MA-0033

Tissue Grossing Table without sink with downdraft extraction. For wall positioning. Mainly used for gross anatomy after dissecting. Versions for 1 or more persons.

Type: „0033“

Description and features of standard version

Most modern tissue grossing concept for guaranteed pollution free (formalin) air in breathing zone of the user. Connection to in house extraction system.

- solid 40mm stainless steel frame construction
- height adjustable feet (40mm).
- 3-part perforated working area (1100Wx600Dmm) w. finger holes.
- grossing board (material: corian) with black „macro field“, 100mm dia., for optimum differentiation). Size: 500 x 400mm Incl. 4 feet in order to assure optimum extraction
- 2 storage platforms below working area for literature and notes
- 1 extraction tube downwards through the ceiling
- front of table is covered (like photo on top)

Optional accessories

- magnifier (3:1) with cold light and flexible support

Technical data

Material: complete stainless steel EN DIN 1.4301 double polished

Material: working surface made of type 316 (4571) stainless steel

Dim.: 1250(W) x 600(D) x 750(H)mm

Height of working area: 750mm (sitting position)

Extraction-air volume: 200-500m³/h adjustable

Air velocity: 0,25m/sec. over work area

Vacuum: 100 Pasc.

Duct: tube 150mm dia.



MA-0033 Tissue Grossing Table with 90° lip for wall protection and downdraft extraction



MA-0033 Tissue Grossing Tables with downdraft extraction for students

Product Specification



Country of Origin

Manufactured by MEDIS MT GmbH in Germany acc. to EN and German norms and the safety regulations of the EC and Germany

